

Microbiology-Molecular Swab Collection Devices for Laboratory Tests 38.15.A03

Please confirm collection device has not expired before using

Test Description	Code	Collection Device	Lawson #
Culture Bacterial, Other with Gram Stain (Sterile container or ESwab)	MICR0032	<u>PREFERRED</u> - Sterile	
Culture Bacterial, Other (Sterile container or ESwab)	MICR0003	screw-cap container with a portion of the specimen (tissue, fluids)	
Culture Bacterial, Anaerobe (Sterile container or ESwab)	MICR0009		
Culture, Screen for Single Pathogen (e.g. VRE) (Sterile container or ESwab)	MICR0005		
Culture, Yeast (Sterile container or ESwab)	MICR5007	Examples of sterile containers	6264320
Culture Bacterial, Respiratory (<u>Nasal and Throat Sources Only</u>) (<u>Sterile container or ESwab)</u>	MICR0001	Alternate option - ESwab (<u>Only</u> if unable to send portion of specimen)	6243783
Culture Bacterial, Respiratory (Nasopharyngeal Sources Only) (ESWAB COLLT FLEX MINITIP ONLY)	Wildings 1	ESWAB REG WHT	
Culture Bacterial, Respiratory with Gram Stain (e.g. Sputum) (Sterile container only)	MICR0031		
Culture, Acid Fast Bacillus (Sterile container only)	MICR0019	ESWAB FLEX MINITIP	
Culture, Fungal (Sterile container only)	MICR0023	6243783 (Sources requiring smaller swab e.g.	
Culture Bacterial, Orthopedic Extended with Gram Stain (Sterile container only)	MICR0046	NP, Sinus)	
Culture, Neisseria (Non-Genital Sources Only)	MICR0044		
Group A Streptococcus by NAT	NBLD0533	★・川子 三回松 平 	6264320
Wet Prep	NBLD0305	KT COLLT COPAN ESWAB REG WHT	
Group B Strep, NAD, No Penicillin Allergy	NBLD5000	1 ** Control (((((((((((((((((((6117712
Group B Strep, NAD Penicillin Allergy, Reflex to Susceptibility	NBLD5011	SWAB CLTR DBL LIQ STUART	
MRSA Nasal Screen, NAT (SWAB CLTR DBL LIQ STUART 6117712 also accepted)	NBLD0394		6264320
Staph aureus/MRSA NAT, Superficial Wound (SWAB CLTR DBL LIQ STUART 6117712 also accepted)	NBLD0494	KT COLLT COPAN ESWAB REG WHT	
B Pertussis/Parapertussis Nucleic Acid Detection	NBLD0274	G. O. D. Helward Vird Transport	
Comprehensive Respiratory Pathogen Panel, NAD	NBLD0481	Mycolgans and Uraplasmas	
Herpes Simplex Virus Nucleic Acid Detection, Non-CSF	NBLD0132	Control Contro	
HSV and VZV Nucleic Acid Detection, Non-CSF	NBLD0360	UTM 3ML FLEX MINI TIP SWAB	6279645
Influenza A and B and RSV Nucleic Acid Detection	NBLD0375	_	6260951
Influenza A/B Nucleic Acid Detection	NBLD0381		
RSV by Nucleic Acid Detection	NBLD0519	MEDIA IIVT W ELOCK FLEX 2M	
Varicella Zoster Virus Nucleic Acid Detection, Non-CSF	NBLD036	MEDIA UVT W FLOCK FLEX 3ML	



Microbiology-Molecular Swab Collection Devices for Laboratory Tests 38.15.A03

Please confirm collection device has not expired before using

Test Description	Code	Collection Device	Lawson #
Chlamydia and Neisseria gonorrhea Nucleic Acid Detection – Vaginal/Rectal/Oropharyngeal	NBLD0267	WINDLE MESSONSTANCE AND THE	
Chlamydia Trachomatis Nucleic Acid Detection – Vaginal/Rectal/Oropharyngeal	NBLD0294	The state of the s	6295646
Neisseria Gonorrhea Nucleic Acid Detection – Vaginal/Rectal/Oropharyngeal	NBLD0295	KT SAMP SWAB PCR MEDIA UNI	
Chlamydia and Neisseria gonorrhea Nucleic Acid Detection - Urine	NBLD0267		
Chlamydia Trachomatis Nucleic Acid Detection - Urine	NBLD0294	Have the second	6246093
Neisseria Gonorrhea Nucleic Acid Detection - Urine	NBLD0295	KT SAMP COBAS PCR URN	
Chlamydia and Neisseria gonorrhea Nucleic Acid Detection – Other (sources other than Vaginal/ Rectal/Oropharyngeal or Urine)	NBLD0615	ACCOUNTY OF THE PARTY OF THE PA	Call the lab
Trichomonas Vaginalis by Nucleic Acid Detection – Vaginal/Endocervical	NBLD0545	XPERT VAG ENDOCERV SWAB	6303032
Trichomonas Vaginalis by Nucleic Acid Detection - Urine	NBLD0545	Xyper to response of chester 80	6272691
Multiplex Vaginal Panel, NAT	NBLD0654	XPERT VAG ENDOCERV SWAB	6303032