




Microbiology-Molecular Swab Collection Devices for Laboratory Tests 38.15.A03

Please confirm collection device has not expired before using

Test Description	Code	Collection Device	Lawson #
Culture Bacterial, Other with Gram Stain (<u>Sterile container or ESwab</u>)	MICR0032	<p style="color: red; text-align: center;">PREFERRED - Sterile screw-cap container with a portion of the specimen (tissue, fluids)</p> <div style="display: flex; justify-content: space-around; align-items: center;">   </div> <p style="color: red; text-align: center;"><i>Examples of sterile containers</i></p> <p>Alternate option - ESwab (Only if unable to send portion of specimen)</p> <div style="display: flex; justify-content: space-around; align-items: center;">  <div> <p>ESWAB REG WHT 6264320</p> </div> </div> <div style="display: flex; justify-content: space-around; align-items: center;">  <div> <p>ESWAB FLEX MINITIP 6243783</p> <p>(Sources requiring smaller swab e.g. NP, Sinus)</p> </div> </div>	6264320 6243783
Culture Bacterial, Other (<u>Sterile container or ESwab</u>)	MICR0003		
Culture Bacterial, Anaerobe (<u>Sterile container or ESwab</u>)	MICR0009		
Culture, Screen for Single Pathogen (e.g. VRE) (<u>Sterile container or ESwab</u>)	MICR0005		
Culture, Yeast (<u>Sterile container or ESwab</u>)	MICR5007		
Culture Bacterial, Upper Respiratory (Nasal and Throat Only) (<u>Sterile container or ESwab</u>)	MICR0001	<p>ESWAB REG WHT 6264320</p> <p>ESWAB FLEX MINITIP 6243783</p> <p>(Sources requiring smaller swab e.g. NP, Sinus)</p>	6264320 6243783
Culture Bacterial, Upper Respiratory (Nasopharyngeal Only) (ESWAB COLLT FLEX MINITIP ONLY)			
Culture Bacterial, Lower Respiratory with Gram Stain (e.g. Sputum) (<u>Sterile container only</u>)	MICR0031		
Culture, Acid Fast Bacillus (<u>Sterile container only</u>)	MICR0019		
Culture, Fungal (<u>Sterile container only</u>)	MICR0023		
Culture Bacterial, Orthopedic Extended with Gram Stain (<u>Sterile container only</u>)	MICR0046	 <p>KT COLLT COPAN ESWAB REG WHT</p>	6264320
Culture, Neisseria (Non-Genital Sources Only)	MICR0044		
Group A Streptococcus by NAT	NBLD0533	 <p>SWAB CLTR DBL LIQ STUART</p>	6117712
Group B Strep by NAD, No Penicillin Allergy	NBLD5000		
Group B Strep NAD Penicillin Allergy, Reflex to Susceptibility	NBLD5011	 <p>KT COLLT COPAN ESWAB REG WHT</p>	6264320
MRSA Nasal Screen, NAT (SWAB CLTR DBL LIQ STUART 6117712 also accepted)	NBLD0394		
MRSA Nasal Screen, NAT Reflex to Screen for Single Pathogen (SWAB CLTR DBL LIQ STUART 6117712/ESWAB FLEX MINITIP 6243783 also accepted)	NBLD0699		
Staph aureus/MRSA NAT, Superficial Wound (SWAB CLTR DBL LIQ STUART 6117712 also accepted)	NBLD0494		




Microbiology-Molecular Swab Collection Devices for Laboratory Tests 38.15.A03

Please confirm collection device has not expired before using

Test Description	Code	Collection Device	Lawson #
B Pertussis/Parapertussis Nucleic Acid Detection	NBLD0274	 UTM 3ML FLEX MINI TIP SWAB	6279645 6260951
Comprehensive Respiratory Pathogen Panel, NAD	NBLD0481		
Herpes Simplex Virus Nucleic Acid Detection, Non-CSF	NBLD0132		
HSV and VZV Nucleic Acid Detection, Non-CSF	NBLD0360		
Influenza A and B and RSV Nucleic Acid Detection	NBLD0375		
Influenza A/B Nucleic Acid Detection	NBLD0381		
RSV by Nucleic Acid Detection	NBLD0519		
SARS-CoV-2, Influenza A+B, and RSV NAT Panel	NBLD0630		
SARS-CoV-2 RNA, Qualitative RT-PCR	NBLD0624		
Varicella Zoster Virus Nucleic Acid Detection, Non-CSF	NBLD036		
Chlamydia and Neisseria gonorrhea Nucleic Acid Detection – Other (Sources other than Vaginal, Rectal, Oropharyngeal, and Urine)	NBLD0615		Call the lab
Chlamydia and Neisseria gonorrhea NAD – Vaginal/Rectal/Oropharyngeal Urgent Care and Outpatients	NBLD0267	 KT SAMP SWAB PCR MEDIA UNI Urgent Care and Outpatients	6295646
Chlamydia Trachomatis Nucleic Acid Detection – Vaginal/Rectal/Oropharyngeal Urgent Care and Outpatients	NBLD0294		
Neisseria Gonorrhea Nucleic Acid Detection – Vaginal/Rectal/Oropharyngeal Urgent Care and Outpatients	NBLD0295		

Microbiology-Molecular Swab Collection Devices for Laboratory Tests 38.15.A03

Please confirm collection device has not expired before using

Test Description	Code	Collection Device	Lawson #
Chlamydia and Neisseria Gonorrhea NAD –Urine Urgent Care and Outpatients	NBLD0267	 KT SAMP COBAS PCR URN Urgent Care and Outpatients	6246093
Chlamydia Trachomatis Nucleic Acid Detection – Urine Urgent Care and Outpatients	NBLD0294		
Neisseria Gonorrhea Nucleic Acid Detection – Urine Urgent Care and Outpatients	NBLD0295		
Chlamydia and Neisseria gonorrhea Nucleic Acid Detection – Vaginal/Endocervical Emergency Department and Inpatients Only	NBLD0545	 XPert SWAB SPECIMEN COLLECTION KIT	6329043
Chlamydia trachomatis Nucleic Acid Detection – Vaginal/Endocervical Emergency Department and Inpatients Only	NBLD0294		
Neisseria gonorrhea Nucleic Acid Detection – Vaginal/Endocervical Emergency Department and Inpatients Only	NBLD0295		
Multiplex Vaginal Panel, NAD	NBLD0654		
Trichomonas Vaginalis by Nucleic Acid Detection – Vaginal/Endocervical	NBLD0545		
Chlamydia and Neisseria gonorrhea Nucleic Acid Detection – Vaginal/Endocervical Emergency Department and Inpatients Only	NBLD0545	 XPert URINE KT COLLT XPert URN	6272691
Chlamydia trachomatis Nucleic Acid Detection – Vaginal/Endocervical Emergency Department and Inpatients Only	NBLD0294		
Neisseria gonorrhea Nucleic Acid Detection – Vaginal/Endocervical Emergency Department and Inpatients Only	NBLD0295		
Trichomonas Vaginalis by Nucleic Acid Detection - Urine	NBLD0545		